

Brochure 2026 Uromat HD



Uromat HD

Key Features

- Easy access from all sides thanks to its small footprint and curved X-ray column
- Isocentric motorized tilting
- Optimal working height for Urologists and Nurses
- 42 x 43 cm Canon B1 detector
- Up to 30 frames per second
- 65 kW or 80 kW high frequency generator
- Canon CXDI RF software and post-processing
- Low dose imaging
- Static and dynamic imaging
- Highly configurable with modular design
- Multifunctional footswitch
- Easy to clean
- Minimal installation space required
- Possibility to move x-ray tube into parking position



Uromat HD

New Detector & Image Quality



Canon CXDI-RF/NE software

- Intuitive to use
- Static and dynamic imaging
- Comprehensive image processing
- Full DICOM package



Canon B1 Detector

- 42 x 43 cm dynamic flat panel detector
- Up to 30 frames per second
- 160 micron pixel size

Uromat HD

Options & Specifications



Uromat HD

Options & Specifications

- **Table Top extension**
- **Foot rest for standing examinations (only in combination with table top extension)**
- **Arm support**
- **Monitor Support Arm for up to two monitors**
- **Additional Support Arm with X-ray protective shield and lamp**
- **Drain Bag Mount**
- **Paper roll holder**
- **Compression Band**
- **Micturition Seat incl Foot Rest**
- **Headrest with Holder**
- **Leg supports**
- **Intercom**



Technical Data	Uromat HD
Table Height	64 - 116,5 cm
Tilting	Up to 88° for standing examinations and up to 20° Trendelenburg
Focus-Detector Distance	115 cm
Table top size	120 x 76 cm
Table top size with table top extension	198 x 76 cm
Longitudinal movement of the table top	± 24 cm
Transverse movement of the table top	± 7 cm
Detector size	42 x 43 cm
High Frequency Generator	65 or 80 kW
DAP	Integrated DAP meter
Anti-scattered grid	Removable
Pulse rate	0,5 - 30 frames per second
Software	Canon CXDI RF
Minimal installation space	370 x 210 x 260 cm (LxWxH)



**Oldelft
Benelux**

A CANON COMPANY

Wiltonstraat 41, 3905 KW Veenendaal, The Netherlands
T +31 318 583 400 E info@delftdi.com I www.oldelftbenelux.nl

Canon
CANON GROUP